

SN-50F6

Syringe pump

► International standard syringe pump

- **Compatible:** Automatic syringe calibration, suitable for all brands
- **Convenient:** "Openable slider hook", convenient to clamp syringe plunger
- **Powerful:** Battery capable to support more than 6 hours for one channel continuous working
- **Energy saving:** Idle channel can be set to sleeping mode



Sino Medical-Device Technology Co., Ltd.

www.sinomdt.com

Product Registration No.: Year 2009, number 2540092, Guangdong FDA

Patent Acceptance No.: 200830155639.2

SN-50F6

Syringe Pump



SPECIFICATION (Dual Channel)

Syringes	10, 20, 30, 50 mL of calibrated syringes	
Flow rate	0.1~1500 mL/h (specified value 0.1 999.9 mL/h, 0.1 mL/h step) (specified value 1000 mL/h or above, 1mL/h step) 0.1~ 900 mL/h (0.1 mL/h step) 0.1~ 600 mL/h (0.1 mL/h step) 0.1~ 400 mL/h (0.1 mL/h step)	50 mL syringe 30 mL syringe 20 mL syringe 10 mL syringe
Purge rate	1500 mL/h 900 mL/h 600 mL/h 400 mL/h	50 mL syringe 30 mL syringe 20 mL syringe 10 mL syringe
Accuracy of flow rate (Nominal value)	Accuracy including syringe: within $\pm 2\%$ Mechanical accuracy: within $\pm 1\%$	
Display range of volume delivered	0.1~ 9999 mL (0.1~999.9 mL, 0.1 mL step) (1000 ~ 9999 mL, 1 mL step)	
Setting of delivery limit	0.1~ 9999 mL (0.1~999.9 mL, 0.1 mL step) (1000 ~ 9999 mL, 1 mL step)	
Occlusion detection pressure	'H': 800 $\pm$ 200 mmHg (106.7 $\pm$ 26.7 kPa) or 1.09 $\pm$ 0.27 kgf/cm ² 'C': 500 $\pm$ 100 mmHg (66.7 $\pm$ 13.3 kPa) or 0.68 $\pm$ 0.13 kgf/cm ² 'L': 300 $\pm$ 100 mmHg (40.7 $\pm$ 13.3 kPa) or 0.41 $\pm$ 0.13 kgf/cm ²	
Alarms	'Occlusion', 'Nearly empty', 'Finish', 'Syringe dislocated', 'Plunger/Clutch disengaged', 'Battery low', 'Battery exhausted', 'Power cable disconnected', 'Flow rate over limit', 'Volume over limit', 'System error', 'Non-operation time over'	
Power	AC110~220V $\pm 10\%$, 50~60 Hz ± 1 Hz Inner NiMH battery: DC12 V Continuous operation time of battery: over 6 hours for one channel, over 4 hours for dual channel when battery is recharged for 12 hours, at flow rate of 5 mL/h.	
Operation conditions	Ambient temperature: 5 to 40 °C, Relative humidity: 20% to 80%	
Storage conditions	Ambient temperature: -40 to 55 °C, Relative humidity: $\leq 95\%$	
Electric shock protection/class	Class I and internal power supply: Type CF, Continuous operation	
Water proof level	IPX4	
Weight	3.45kg (with clamp base)	
Dimension	354mm (Width) $\times$ 190mm (Height) $\times$ 128mm (Depth)	

Function

Dual-channel with different flow rate	Operation error prompt
Channel sleeping function	Fast rate control
Identify syringe size automatically	"Purge" button activate safe protection
3-level occlusion alarm limits setting	Flow rate setting button lock
Limits setting	Abundant optic and sound alarm
Automatically reduce impact dose release after occlusion	AC and internal power indication
Delivered volume inquiry	Standard RS232 interface
KVO flow rate	500 history records

Manufacturer: Sino Medical-Device Technology Co., Ltd.

Add: Room 2-D-2, Huake Building, Oriental Hi-Tech Industrial Park,
OCT Nanshan District, Shenzhen, P.R.China Zip: 518053

Tel: (86)755-26755692 Fax: (86)755-26755687

Website: <http://www.sinomdt.com>

Email: information@sinomdt.com

Agent: Hangzhou Santtel Electric Technologies Co., Ltd.

Add: Room A405 Technology Building, No. 90 Wensan Road,
Hangzhou, P.R.China Zip: 310012

Tel: (86)571-81951038 Fax: (86)571-81951037

Website: <http://www.sinomdt.com>

Email: information@sinomdt.com

This document is for reference only, please refer to product user manual for detail, products update without notice.